



CLEARWINGS® BRACKET

MBT Prescription .018

Tooth	Torque	Angulation	Wide	Right	Left
Upper - Red					
				Item #	Item #
Central	17°	5°	3.20	406-11	406-21
Lateral	10°	9°	2.95	406-12	406-22
Cuspid	0°	8°	3.20	406-13	406-23
Cuspid - Hook	0°	8°	3.20	406-13H	406-23/H
1 Bicuspid	-7°	0°	3.20	406-14	406-24
1 Bicuspid - hook	-7°	0°	3.20	406-14H	406-24H
2 Bicuspid	-7°	0°	3.20	406-15	406-15
2 Bicuspid - hook	-7°	0°	3.20	406-15H	406-15H
Lower - Black					
Lower Anteriors	-6°	0°	2.95	406-41/42	406-31/32
Cuspid	0°	3°	3.50	406-43	406-33
Cuspid - Hook	-0°	3°	3.20	406-43/H	406-33/H
1 Bicuspid	-12°	2°	3.20	406-44	406-34
1 Bicuspid - hook	-12°	2°	3.20	406-44H	406-34H
2 Bicuspid	-17°	2°	3.20	406-45	406-45
2 Bicuspid - hook	-17°	2°	3.20	406-45H	406-45H

MBT Prescription .022

Tooth	Torque	Angulation	Wide	Right	Left
Upper - Red					
				Item #	Item #
Central	17°	5°	3.20	466-11	466-21
Lateral	10°	9°	2.95	466-12	466-22
Cuspid	0°	8°	3.20	466-13	466-23
Cuspid - Hook	0°	8°	3.20	466-13H	466-23/H
1 Bicuspid	-2°	0°	3.20	466-14	466-24
1 Bicuspid - hook	-2°	0°	3.20	466-14H	466-24H
2 Bicuspid	-7°	0°	3.20	466-15	466-15
2 Bicuspid - hook	-7°	0°	3.20	466-15H	466-15H
Lower - Black					
Lower Anteriors	-6°	0°	2.95	466-41/42	466-31/32
Cuspid	0°	3°	3.50	466-43	466-33
Cuspid - Hook	-0°	3°	3.20	466-43/H	466-33/H
1 Bicuspid	-12°	2°	3.20	466-44	466-34
1 Bicuspid - hook	-12°	2°	3.20	466-44H	466-34H
2 Bicuspid	-17°	2°	3.20	466-45	466-45
2 Bicuspid - hook	-17°	2°	3.20	466-45H	466-45H



CLEARWINGS® BRACKET

Roth Prescription .018

Tooth	Torque	Angulation	Wide	Right	Left
Upper - Red					
				Item #	Item #
Central	11°	5°	3.20	305-11	305-21
Lateral	8°	9°	2.95	305-12	305-22
Cuspid	-2°	8°	3.20	305-13	305-23
Cuspid - Hook	-2°	8°	3.20	305-13H	305-23/H
1 Bicuspid	-7°	0°	3.20	305-14	305-24
1 Bicuspid - hook	-7°	0°	3.20	305-14H	305-24H
2 Bicuspid	-7°	0°	3.20	305-15	305-15
2 Bicuspid - hook	-7°	0°	3.20	305-15H	305-15H
Lower - Black					
Lower Anteriors	0°	0°	2.95	305-41/42	305-31/32
Cuspid	-11°	7°	3.50	305-43	305-33
Cuspid - Hook	-11°	7°	3.20	305-43/H	305-33/H
1 Bicuspid	-17°	3°	3.20	305-44	305-34
1 Bicuspid - hook	-17°	3°	3.20	305-44H	305-34H
2 Bicuspid	-21°	6°	3.20	305-45	305-45
2 Bicuspid - hook	-21°	6°	3.20	305-45H	305-45H

Roth Prescription .022

Tooth	Torque	Angulation	Wide	Right	Left
Upper - Red					
				Item #	Item #
Central	11°	5°	3.20	355-11	355-21
Lateral	8°	9°	2.95	355-12	355-22
Cuspid	-2°	8°	3.20	355-13	355-23
Cuspid - Hook	-2°	8°	3.20	355-13H	355-23/H
1 Bicuspid	-7°	0°	3.20	355-14	355-24
1 Bicuspid - hook	-7°	0°	3.20	355-14H	355-24H
2 Bicuspid	-7°	0°	3.20	355-15	355-15
2 Bicuspid - hook	-7°	0°	3.20	355-15H	355-15H
Lower - Black					
Lower Anteriors	0°	0°	2.95	355-41/42	355-31/32
Cuspid	-11°	7°	3.50	355-43	355-33
Cuspid - Hook	-11°	7°	3.20	355-43/H	355-33/H
1 Bicuspid	-17°	3°	3.20	355-44	355-34
1 Bicuspid - hook	-17°	3°	3.20	355-44H	355-34H
2 Bicuspid	-21°	6°	3.20	355-45	355-45
2 Bicuspid - hook	-21°	6°	3.20	355-45H	355-45H