

CROWN[®] BRACKET



Roth Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red					Item #	Item #
Central	12°	5°	0.84	2.85	105-11	105-21
Lateral	8°	9°	1.15	2.55	105-12	105-22
Cuspid	-2°	9°	0.60	2.90	105-13	105-23
Cuspid - Hook	-2°	9°	0.60	2.90	105-13H	105-23/H
1-2 Bicupids	-7°	0°	0.68	2.95	105-14/15	105-24/25
1-2 Bicuspids - hook	-7°	0°	0.68	2.95	105-14/15/H	105-24/25/H
Lower - Black						
Lower Anteriors	-1°	0°	1.05	2.65	105-41/42	105-31/32
Cuspid	-1°	7°	0.56	2.90	105-43	105-33
Cuspid - Hook	-11°	7°	0.56	2.90	105-43/H	105-33/H
1 Bicuspid	-17°	0°	0.54	2.95	105-44	105-34
2 Bicuspid	-22°	0°	0.52	2.95	105-45	105-35
1 Bicuspid - hook	-17°	0°	0.54	2.95	105-44/H	105-34/H
2 Bicuspid - hook	-22°	0°	0.52	2.95	105-45/H	105-35/H

Roth Prescription .022

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red					Item #	Item #
Central	7°	5°	0.79	3.45	155-11	155-21
Lateral	3°	9°	1.28	2.90	155-12	155-22
Cuspid	-7°	11°	0.69	3.15	155-13	155-23
Cuspid - Hook	-7°	11°	0.69	3.15	155-13H	155-23/H
1-2 Bicupids	-7°	2°	0.69	3.20	155-14/15	155-24/25
1-2 Bicuspids - hook	-7°	2°	0.69	3.20	155-14/15/H	155-24/25/H
Lower - Black						
Lower Anteriors	-1°	2°	1.28	2.35	155-41/42	155-31/32
Cuspid	-11°	5°	0.64	3.15	155-43	155-33
Cuspid - Hook	-11°	5°	0.64	3.15	155-43/H	155-33/H
1 Bicuspid	-17°	2°	0.46	3.20	155-44	155-34
2 Bicuspid	-22°	2°	0.46	3.20	155-45	155-35
1 Bicuspid - hook	-17°	2°	0.46	3.20	155-44/H	155-34/H
2 Bicuspid - hook	-22°	2°	0.46	3.20	155-45/H	155-35/H