




ECONOLINE® BRACKET

Roth Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
						Item #
Upper - Red						
Central	7°	5°	0.79	3.45	05-11	05-21
Lateral	3°	9°	1.28	2.90	05-12	05-22
Cuspid	-7°	11°	0.69	3.15	05-13	05-23
Cuspid - Hook	-7°	11°	0.69	3.15	05-13H	05-23/H
1-2 Bicuspids	-7°	2°	0.69	3.20	05-14/15	05-24/25
1-2 Bicuspid - hook	-7°	2°	0.69	3.20	05-14/15/H	05-24/25/H
Lower - Black						
Lower Anteriors	-1°	2°	1.28	2.35	05-41/42	05-31/32
Cuspid	-11°	5°	0.64	3.15	05-43	05-33
Cuspid - Hook	-11°	5°	0.64	3.15	05-43/H	05-33/H
1 Bicuspid	-17°	2°	0.46	3.20	05-44	05-34
2 Bicuspid	-22°	2°	0.46	3.20	05-45	05-35
1 Bicuspid - hook	-17°	2°	0.46	3.20	05-44/H	05-34/H
2 Bicuspid - hook	-22°	2°	0.46	3.20	05-45/H	05-35/H

Roth Prescription .022

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
						Item #
Upper - Red						
Central	7°	5°	0.79	3.45	55-11	55-21
Lateral	3°	9°	1.28	2.90	55-12	55-22
Cuspid	-7°	11°	0.69	3.15	55-13	55-23
Cuspid - Hook	-7°	11°	0.69	3.15	55-13H	55-23/H
1-2 Bicuspids	-7°	2°	0.69	3.20	55-14/15	55-24/25
1-2 Bicuspid - hook	-7°	2°	0.69	3.20	55-14/15/H	55-24/25/H
Lower - Black						
Lower Anteriors	-1°	2°	1.28	2.35	55-41/42	55-31/32
Cuspid	-11°	5°	0.64	3.15	55-43	55-33
Cuspid - Hook	-11°	5°	0.64	3.15	55-43/H	55-33/H
1 Bicuspid	-17°	2°	0.46	3.20	55-44	55-34
2 Bicuspid	-22°	2°	0.46	3.20	55-45	55-35
1 Bicuspid - hook	-17°	2°	0.46	3.20	55-44/H	55-34/H
2 Bicuspid - hook	-22°	2°	0.46	3.20	55-45/H	55-35/H




ECONOLINE® BRACKET

Andrews Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
						Item #
Upper - Red						
Central	7°	5°	0.79	3.45	03-11	03-21
Lateral	3°	9°	1.28	2.90	03-12	03-22
Cuspid	-7°	11°	0.69	3.15	03-13	03-23
Cuspid - Hook	-7°	11°	0.69	3.15	03-13H	03-23/H
1-2 Bicuspids	-7°	2°	0.69	3.20	03-14/15	03-24/25
1-2 Bicuspid - hook	-7°	2°	0.69	3.20	03-14/15/H	03-24/25/H
Lower - Black						
Lower Anteriors	-1°	2°	1.28	2.35	03-41/42	03-31/32
Cuspid	-11°	5°	0.64	3.15	03-43	03-33
Cuspid - Hook	-11°	5°	0.64	3.15	03-43/H	03-33/H
1 Bicuspid	-17°	2°	0.46	3.20	03-44	03-34
2 Bicuspid	-22°	2°	0.46	3.20	03-45	03-35
1 Bicuspid - hook	-17°	2°	0.46	3.20	03-44/H	03-34/H
2 Bicuspid - hook	-22°	2°	0.46	3.20	03-45/H	03-35/H

Andrews Prescription .022

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
						Item #
Upper - Red						
Central	7°	5°	0.79	3.45	33-11	33-21
Lateral	3°	9°	1.28	2.90	33-12	33-22
Cuspid	-7°	11°	0.69	3.15	33-13	33-23
Cuspid - Hook	-7°	11°	0.69	3.15	33-13H	33-23/H
1-2 Bicuspids	-7°	2°	0.69	3.20	33-14/15	33-24/25
1-2 Bicuspid - hook	-7°	2°	0.69	3.20	33-14/15/H	33-24/25/H
Lower - Black						
Lower Anteriors	-1°	2°	1.28	2.35	33-41/42	33-31/32
Cuspid	-11°	5°	0.64	3.15	33-43	33-33
Cuspid - Hook	-11°	5°	0.64	3.15	33-43/H	33-33/H
1 Bicuspid	-17°	2°	0.46	3.20	33-44	33-34
2 Bicuspid	-22°	2°	0.46	3.20	33-45	33-35
1 Bicuspid - hook	-17°	2°	0.46	3.20	33-44/H	33-34/H
2 Bicuspid - hook	-22°	2°	0.46	3.20	33-45/H	33-35/H




ECONOLINE® BRACKET

MBT Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red						
					Item #	Item #
Central	17°	4°	0.79	3.45	06-11	06-21
Lateral	10°	8°	2.90	2.90	06-12	06-22
Cuspid	-7°	8°	3.15	3.15	06-13	06-23
Cuspid - Hook	0°	8°	3.15	3.15	06-13H	06-23/H
1 Bicuspid	-7°	0°	3.20	3.20	06-14	06-24
2 Bicuspid	-7°	0°	3.20	3.20	06-15	06-25
1 Bicuspid - hook	-7°	0°	3.20	3.20	06-14/H	06-24/H
1 Bicuspid - hook	-7°	0°	3.20	3.20	06-15/H	06-25/H
Lower - Black						
Lower Anteriors	-6°	0°	1.20	2.35	06-41/42	06-31/32
Cuspid	-6°	3°	0.80	3.15	06-43	06-33
Cuspid - Hook	-0°	5°	0.80	3.15	06-43/H	06-33/H
1 Bicuspid	-12°	2°	0.80	3.20	06-44	06-34
2 Bicuspid	-17°	2°	0.80	3.20	06-45	06-35
1 Bicuspid - hook	-12°	2°	0.80	3.20	06-44/H	06-34/H
2 Bicuspid - hook	-17°	2°	0.80	3.20	05-45/H	06-35/H

MBT Prescription .022

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red						
					Item #	Item #
Central	17°	4°	0.79	3.45	66-11	66-21
Lateral	10°	8°	2.90	2.90	66-12	66-22
Cuspid	-7°	8°	3.15	3.15	66-13	66-23
Cuspid - Hook	0°	8°	3.15	3.15	66-13H	66-23/H
1 Bicuspid	-7°	0°	3.20	3.20	66-14	66-24
2 Bicuspid	-7°	0°	3.20	3.20	66-15	66-25
1 Bicuspid - hook	-7°	0°	3.20	3.20	66-14/H	66-24/H
1 Bicuspid - hook	-7°	0°	3.20	3.20	66-15/H	66-25/H
Lower - Black						
Lower Anteriors	-6°	0°	1.20	2.35	66-41/42	66-31/32
Cuspid	-6°	3°	0.80	3.15	66-43	66-33
Cuspid - Hook	-0°	5°	0.80	3.15	66-43/H	66-33/H
1 Bicuspid	-12°	2°	0.80	3.20	66-44	66-34
2 Bicuspid	-17°	2°	0.80	3.20	66-45	66-35
1 Bicuspid - hook	-12°	2°	0.80	3.20	66-44/H	66-34/H
2 Bicuspid - hook	-17°	2°	0.80	3.20	65-45/H	66-35/H



ECONOLINE® BRACKET

Standard Edgewise Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right & Left
Upper - Red					
					Item #
Central	0°	0°	0.75	3.45	01-11/21
Lateral	0°	0°	0.75	2.90	01-12/22
Cuspid U - L	0°	0°	0.75	3.15	01-13/23/33/43/23/33/4313/23/33/43
Cuspid - Hook UR - LL	0°	0°	0.75	3.15	01-13/33/H
Cuspid - Hook UL - LR	0°	0°	0.75	3.15	01-23/43/H
1-2 Bicuspids U - L	0°	0°	0.75	3.20	01-14/24/34/44
1-2 Bicuspids - hook R	0°	0°	0.75	3.20	01-14/15/H
1-2 Bicuspids - hook L	0°	0°	0.75	3.20	01-24/25/H
Lower - Black					
Lower Anterior	0°	0°	0.75	2.35	01-31/32/41/42
1-2 Bicuspids - hook R	0°	0°	0.75	3.20	01-34/35/H
1-2 Bicuspids - hook L	0°	0°	0.75	3.20	01-44/45/H

Standard Edgewise Prescription .022

Tooth	Torque	Angulation	In/Out	Wide	Right & Left
Upper - Red					
					Item #
Central	0°	0°	0.75	3.45	11-11/21
Lateral	0°	0°	0.75	2.90	11-12/22
Cuspid U - L	0°	0°	0.75	3.15	11-13/23/33/43/23/33/4313/23/33/43
Cuspid - Hook UR - LL	0°	0°	0.75	3.15	11-13/33/H
Cuspid - Hook UL - LR	0°	0°	0.75	3.15	11-23/43/H
1-2 Bicuspids U - L	0°	0°	0.75	3.20	11-14/24/34/44
1-2 Bicuspids - hook R	0°	0°	0.75	3.20	11-14/15/H
1-2 Bicuspids - hook L	0°	0°	0.75	3.20	11-24/25/H
Lower - Black					
Lower Anterior	0°	0°	0.75	2.35	11-31/32/41/42
1-2 Bicuspids - hook R	0°	0°	0.75	3.20	11-34/35/H
1-2 Bicuspids - hook L	0°	0°	0.75	3.20	11-44/45/H

ECONOLINE® BRACKET




Ricketts Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red					Item #	Item #
Central	22°	0°	0.80	3.45	02-11	02-21
Lateral	14°	8°	0.80	2.90	02-12	02-22
Cuspid	7°	5°	0.80	3.15	02-13	02-23
Cuspid - Hook	7°	5°	0.80	3.15	02-13H	02-23/H
1 - 2 Bicuspid	0°	0°	0.80	3.20	02-14/15	02-24/25
1 - 2 Bicuspid - hook	0°	0°	0.80	3.20	02-14/15H	02-24/25H
Lower - Black						
Lower Anteriors	0°	0°	0.80	2.35	02-41/42	02-31/32
Cuspid	7°	3°	0.80	3.15	02-43	02-33
Cuspid - Hook	7°	5°	0.80	3.15	02-43/H	02-33/H
1-2 Bicuspid	0°	0°	0.80	3.20	02-44	02-34/35
1-2 Bicuspid - hook	0°	0°	0.80	3.20	02-45	02-34/35/H

Hilgers Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red					Item #	Item #
Central	22°	5°	0.70	3.45	07-11	07-21
Lateral	14°	8°	0.70	2.90	07-12	07-22
Cuspid	7°	7°	0.70	3.15	07-13	07-23
Cuspid - Hook	7°	7°	0.70	3.15	07-13H	07-23/H
1 - 2 Bicuspid	-7°	0°	0.70	3.20	07-14/15	07-24/25
1 - 2 Bicuspid - hook	-7°	0°	0.70	3.20	07-14/15H	07-24/25H
Lower - Black						
Lower Anteriors	-1°	0°	1.28	2.35	07-41/42	07-31/32
Cuspid	7°	5°	0.70	3.15	07-43	07-33
Cuspid - Hook	7°	5°	0.70	3.15	07-43/H	07-33/H
1-Bicuspid	-11°	0°	0.70	3.20	07-44	07-34
2- Bicuspid	-22°	0°	0.70	3.20	07-45	07-35
1- Bicuspid - hook	-11°	0°	0.70	3.20	07-44/H	07-34/H
2- Bicuspid - hook	-22°	0°	0.70	3.20	07-45/H	07-35/H

ECONOLINE® BRACKET



Hasund New 2002 Prescription .018

Tooth	Torque	Angulation	In/Out	Wide	Right	Left
Upper - Red					Item #	Item #
Central	20°	5°	1.10	3.45	04-11	04-21
Lateral	14°	9°	1.28	2.90	04-12	04-22
Cuspid	-2°	9°	0.70	3.15	04-13	04-23
Cuspid - Hook	-2°	9°	0.70	3.15	04-13H	04-23/H
1 - 2 Bicuspid	-10°	2°	0.70	3.20	04-14/15	04-24/25
Lower - Black						
Lower Anteriors	0°	3°	1.10	2.35	04-41/42	04-31/32
Cuspid	7°	0°	0.70	3.15	04-43	04-33
Cuspid - Hook	7°	0°	0.70	3.15	04-43/H	04-33/H
1-Bicuspid	-10°	2°	0.50	3.20	04-44	04-34/35
2- Bicuspid	-15°	2°	0.50	3.20	04-45	04-35

Tweed Technique Prescription .022

Tooth Lower	Torque	Angulation	In/Out	Wide	Item #
Central Upper Right + Left	0°	0°	0.75	2.35	00-1/2
Lateral Upper Right + Left	0°	0°	0.75	2.35	00-1/2
Cuspid Upper + Lower Right + Left	0°	0°	0.75	2.35	003/4/5
1-2 Bicuspid Upper + Lower	0°	0°	0.75	2.35	00/3/4/5
Lower - Black					
1-2 Lower Anteriors	0°	0°	0.75	2.35	00-31/42