

FLAIR® SELF-LIGATING BRACKET



Roth Prescription .018

Tooth Upper	Torque	Angulation	In/Out	Wide	Right	Left
Central	+12°	+5°	0.59	2.6	205M-11	205M-21
Lateral	+8°	+9°	0.94	2.6	205M-12	205M-22
Cuspid	-2°	+9°	0.50	2.8	205M-13	205M-23
1 + 2 Bicuspid	-7°	0°	0.59	2.6	205M-14/25	205M-14/25
Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
1+2 Anteriors	-1°	0°	0.85	2.6	205M-31/42	205M-31/42
Cuspid	-11°	+7°	0.43	2.8	205M-43	205M-33
1 Bicuspid	-17°	0°	0.45	2.8	205M-44	205M-34
2 Bicuspid	-22°	0°	0.47	2.8	205M-45	205M-35

Roth Prescription .022

Tooth Upper	Torque	Angulation	In/Out	Wide	Right	Left
Central	+12°	+5°	0.59	2.6	255M-11	255M-21
Lateral	+8°	+9°	0.94	2.6	255M-12	255M-22
Cuspid	-2°	+9°	0.50	2.8	255M-13	255M-23
1 + 2 Bicuspid	-7°	0°	0.59	2.6	255M-14/25	255M-14/25
Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
1+2 Anteriors	-1°	0°	0.85	2.6	255M-31/42	255M-31/42
Cuspid	-11°	+7°	0.43	2.8	255M-43	255M-33
1 Bicuspid	-17°	0°	0.45	2.8	255M-44	255M-34
2 Bicuspid	-22°	0°	0.47	2.8	255M-45	255M-35

Roth **High Torque** Prescription .022 (more High Torques details on page 20)

Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
Central						
Lateral						
Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
1+2 Anteriors	-6°	0°	0.85	2.6	255M-31/42-6	255M-31/42-6

Roth **High Torque** Prescription .018 (more High Torques details on page 20)

Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
Central						
Lateral						
Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
1+2 Anteriors	-6°	0°	0.85	2.6	205M-31/42-6	205M-31/42-6

FLAIR® SELF-LIGATING BRACKET



MBT Prescription .018

Tooth Upper	Torque	Angulation	In/Out	Wide	Right	Left
Central	+17°	+4°	0.59	2.6	206M-11	206M-21
Lateral	+10°	+8°	0.94	2.6	206M-12	206M-22
Cuspid	-7°	+8°	0.50	2.8	206M-13	206M-23
1 Bicuspid	-7°	0°	0.72	2.8	206M-14/24	206M-14/24
2 Bicuspid	-7°	0°	1.20	2.8	206M-15/25	206M-15/25
Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
1+2 Anteriors	-6°	0°	0.85	2.6	206M-31/42	206M-31/42
Cuspid	-6°	+3°	0.43	2.8	206M-43	206M-33
1 Bicuspid	-12°	+2°	0.45	2.8	206M-44	206M-34
2 Bicuspid	-17°	+2°	0.47	2.8	206M-45	206M-35

MBT Prescription .022

Tooth Upper	Torque	Angulation	In/Out	Wide	Right	Left
Central	+17°	+4°	0.59	2.6	266M-11	266M-21
Lateral	+10°	+8°	0.94	2.6	266M-12	266M-22
Cuspid	-7°	+8°	0.50	2.8	266M-13	266M-23
1 Bicuspid	-7°	0°	0.72	2.8	266M-14/24	266M-14/24
2 Bicuspid	-7°	0°	1.20	2.8	266M-15/25	266M-15/25
Tooth Lower	Torque	Angulation	In/Out	Wide	Right	Left
1+2 Anteriors	-6°	0°	0.85	2.6	266M-31/42	266M-31/42
Cuspid	-6°	+3°	0.43	2.8	266M-43	266M-33
1 Bicuspid	-12°	+2°	0.45	2.8	266M-44	266M-34
2 Bicuspid	-17°	+2°	0.47	2.8	266M-45	266M-35