



NICKELFREE® BRACKET

Roth Prescription .018


Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red					Item #	Item #
Central	+12°	+5°	0.029	0.071	NF05-11	NF05-21
Lateral	+8°	+8°	0.024	0.078	NF05-12	NF05-22
Cuspid	+2°	+13°	0.030	0.065	NF05-13	NF05-23
Cuspid - Hook	+2°	+13°	0.030	0.065	NF05-13H	NF05-23/H
1-2 Bicuspids	-7°	0°	0.031	0.069	NF05-14/15	NF05-24/25
1-2 Bicuspids - hook	-7°	0°	0.031	0.069	NF05-14/15/H	NF05-24/25/H
Lower - Black						
Lower Anteriors	+1°	0°	0.031	0.085	NF05-41/42	NF05-31/32
Cuspid	-11°	+7°	0.029	0.066	NF05-43	NF05-33
Cuspid - Hook	-11°	+7°	0.029	0.066	NF05-43/H	NF05-33/H
1 Bicuspid	-17°	0°	0.029	0.065	NF05-44	NF05-34
2 Bicuspid	-22°	0°	0.029	0.065	NF05-45	NF05-35
1 Bicuspid - hook	-17°	0°	0.029	0.065	NF05-44/H	NF05-34/H
2 Bicuspid - hook	-22°	0°	0.029	0.065	NF05-45/H	NF05-35/H

Roth Prescription .022

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red					Item #	Item #
Central	+12°	+5°	0.029	0.071	NF55-11	NF55-21
Lateral	+8°	+8°	0.024	0.078	NF55-12	NF55-22
Cuspid	+2°	+13°	0.030	0.065	NF55-13	NF55-23
Cuspid - Hook	+2°	+13°	0.030	0.065	NF55-13H	NF55-23/H
1-2 Bicuspids	-7°	0°	0.031	0.069	NF55-14/15	NF55-24/25
1-2 Bicuspids - hook	-7°	0°	0.031	0.069	NF05-14/15/H	NF55-24/25/H
Lower - Black						
Lower Anteriors	+1°	0°	0.031	0.085	NF55-41/42	NF55-31/32
Cuspid	-11°	+7°	0.029	0.066	NF55-43	NF55-33
Cuspid - Hook	-11°	+7°	0.029	0.066	NF55-43/H	NT55-33/H
1 Bicuspid	-17°	0°	0.029	0.065	NF55-44	NF55-34
2 Bicuspid	-22°	0°	0.029	0.065	NF55-45	NF55-35
1 Bicuspid - hook	-17°	0°	0.029	0.065	NF55-44/H	NF55-34/H
2 Bicuspid - hook	-22°	0°	0.029	0.065	NF55-45/H	NF55-35/H

Prescriptions

NICKELFREE® BRACKET



MBT Prescription .018

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red						
					Item #	Item #
Central	+17°	+4°	0.030	0.073	NF06-11	NF06-21
Lateral	+8°	+9°	0.024	0.078	NF06-12	NF06-22
Cuspid	-7°	+8°	0.031	0.063	NF06-13	NF06-23
Cuspid - Hook	-7°	+8°	0.031	0.063	NF06-13H	NF06-23/H
1-2 Bicuspids	-7°	0°	0.031	0.069	NF06-14/15	NF06-24/25
1-2 Bicuspid - hook	-7°	0°	0.031	0.069	NF06-14/15/H	NF06-24/25/H
Lower - Black						
Lower Anteriors	-6°	0°	0.031	0.085	NF06-41/42	NF06-31/32
Cuspid	-6°	+3°	0.029	0.063	NF06-43	NF06-33
Cuspid - Hook	-6°	+3°	0.029	0.063	NF06-43/H	NF06-33/H
1 Bicuspid	-12°	+2°	0.029	0.070	NF06-44	NF06-34
2 Bicuspid	-17°	+2°	0.029	0.070	NF06-45	NF06-35
1 Bicuspid - hook	-12°	+2°	0.029	0.070	NF06-44/H	NF06-34/H
2 Bicuspid - hook	-17°	+2°	0.029	0.070	NF06-45/H	NF06-35/H

MBT Prescription .022

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red						
					Item #	Item #
Central	+17°	+4°	0.030	0.073	NF66-11	NF66-21
Lateral	+8°	+9°	0.024	0.078	NF66-12	NF66-22
Cuspid	-7°	+8°	0.031	0.063	NF66-13	NF66-23
Cuspid - Hook	-7°	+8°	0.031	0.063	NF66-13H	NF66-23/H
1-2 Bicuspids	-7°	0°	0.031	0.069	NF66-14/15	NF66-24/25
1-2 Bicuspid - hook	-7°	0°	0.031	0.069	NF66-14/15/H	NF66-24/25/H
Lower - Black						
Lower Anteriors	-6°	0°	0.031	0.085	NF66-41/42	NF66-31/32
Cuspid	-6°	+3°	0.029	0.063	NF66-43	NF66-33
Cuspid - Hook	-6°	+3°	0.029	0.063	NF66-43/H	NF66-33/H
1 Bicuspid	-12°	+2°	0.029	0.070	NF66-44	NF66-34
2 Bicuspid	-17°	+2°	0.029	0.070	NF66-45	NF66-35
1 Bicuspid - hook	-12°	+2°	0.029	0.070	NF66-44/H	NF66-34/H
2 Bicuspid - hook	-17°	+2°	0.029	0.070	NF66-45/H	NF66-35/H

Prescriptions

NICKELFREE® BRACKET



Hilgers & Ricketts Prescription .018

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red						
					Item #	Item #
Central	+22°	+5°	0.024	0.066	NF02-11	NF02-21
Lateral	+14°	+8°	0.026	0.077	NF02-12	NF02-22
Cuspid	+7°	+7°	0.029	0.061	NF02-13	NF02-23
Cuspid - Hook	+7°	+7°	0.028	0.061	NF02-13H	NF02-23/H
1-2 Bicuspids	-7°	0°	0.031	0.068	NF02-14/15	NF02-24/25
1-2 Bicuspid - hook	-7°	0°	0.031	0.068	NF02-14/15/H	NF02-24/25/H
Lower - Black						
Lower Anteriors	0°	0°	0.031	0.085	NF02-41/42	NF02-31/32
Cuspid	+7°	+5°	0.029	0.058	NF02-43	NF02-33
Cuspid - Hook	+7°	+5°	0.029	0.058	NF02-43/H	NF02-33/H
1 Bicuspid	-11°	0°	0.024	0.066	NF02-44	NF02-34
2 Bicuspid	-17°	0°	0.029	0.065	NF02-45	NV02-35
1 Bicuspid - hook	-11°	0°	0.024	0.066	NF02-44/H	NF02-34/H
2 Bicuspid - hook	-17°	0°	0.029	0.065	NF02-45/H	NF02-35/H

Hilgers & Ricketts Prescription .022

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red						
					Item #	Item #
Central	+22°	+5°	0.024	0.066	NF22-11	NF22-21
Lateral	+14°	+8°	0.026	0.077	NF22-12	NF22-22
Cuspid	+7°	+7°	0.029	0.061	NF22-13	NF22-23
Cuspid - Hook	+7°	+7°	0.028	0.061	NF22-13H	NF22-23/H
1-2 Bicuspids	-7°	0°	0.031	0.068	NF22-14/15	NF22-24/25
1-2 Bicuspid - hook	-7°	0°	0.031	0.068	NF22-14/15/H	NF22-24/25/H
Lower - Black						
Lower Anteriors	0°	0°	0.031	0.085	NF22-41/42	NF22-31/32
Cuspid	+7°	+5°	0.029	0.058	NF22-43	NF22-33
Cuspid - Hook	+7°	+5°	0.029	0.058	NF22-43/H	NF22-33/H
1 Bicuspid	-11°	0°	0.024	0.066	NF22-44	NF22-34
2 Bicuspid	-17°	0°	0.029	0.065	NF22-45	NV22-35
1 Bicuspid - hook	-11°	0°	0.024	0.066	NF22-44/H	NF22-34/H
2 Bicuspid - hook	-17°	0°	0.029	0.065	NF22-45/H	NF22-35/H

NICKELFREE[®] BRACKET

Standard Edgewise Prescription .018

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red					Item #	Item #
Central	0°	0°	0.027	0.069	NF01-11	NF01-21
Lateral	0°	0°	0.024	0.064	NF01-12	NF01-22
Cuspid	0°	0°	0.027	0.056	NF01-13	NF01-23
1-2 Bicupids	0°	0°	0.027	0.056	NF01-14/15	NF01-24/25
Lower - Black						
Lower Anteriors	0°	0°	0.029	0.065	NF01-41/42	NF01-31/32
Cuspid	0°	0°	0.027	0.056	NF01-43	NF01-33
1 Bicuspid	0°	0°	0.027	0.056	NF01-44	NF01-34
2 Bicuspid	0°	0°	0.027	0.056	NF01-45	NF01-35

Standard Edgewise Prescription .022

Tooth	Torque	Angulation	Slot Depth	Height	Right	Left
Upper - Red					Item #	Item #
Central	0°	0°	0.027	0.069	NF11-11	NF11-21
Lateral	0°	0°	0.024	0.064	NF11-12	NF11-22
Cuspid	0°	0°	0.027	0.056	NF11-13	NF11-23
1-2 Bicupids	0°	0°	0.027	0.056	NF11-14/15	NF11-24/25
Lower - Black						
Lower Anteriors	0°	0°	0.029	0.065	NF11-41/42	NF11-31/32
Cuspid	0°	0°	0.027	0.056	NF11-43	NF11-33
1 Bicuspid	0°	0°	0.027	0.056	NF11-44	NF11-34
2 Bicuspid	0°	0°	0.027	0.056	NF11-45	NF11-35