

## NICKEL FREE<sup>®</sup> BUCCAL TUBE



Single Buccal Tube Roth Prescription .018

Tooth	Torque	Offset	In/Out	Right	Left
<b>Upper</b>				Item #	Item #
Single 1st Molar	-10°	+8°	.042	018S10-112R	018S10-112L
Single	0°	+8°	.042	018S10-104R	018S10-104L
Single 1st/2nd Molar	-14°	+8°	.042	018S10-131R	018S10-131L
Universal Single 1st Molar	0°	0°	.028	018S10-001A	018S10-001B
Universal Single	-10°	0°	.028	018S10-011A	018S10-011B
<b>Lower</b>					
Single 1st Molar	-25°	+5°	.042	018S10-252R	018S10-252L
Single 2nd Molar	-20°	0°	.042	018S10-240R	018S10-240L
Single	0°	+5°	.042	018S10-203R	018S10-203L
Universal Single 1st Molar	0°	0°	.028	018S10-001B	018S10-001A
Universal Single	-10°	0°	.028	018S10-011B	018S10-011A



Single Buccal Tube Roth Prescription .022

Tooth	Torque	Offset	In/Out	Right	Left
<b>Upper</b>				Item #	Item #
Single 1st Molar	-10°	+8°	.042	022S10-112R	022S10-112L
Single	0°	+8°	.042	022S10-104R	022S10-104L
Single 1st/2nd Molar	-14°	+8°	.042	022S10-131R	022S10-131L
Universal Single 1st Molar	0°	0°	.028	022S10-001A	022S10-001B
Universal Single	-10°	0°	.028	022S10-011A	022S10-011B
<b>Lower</b>					
Single 1st Molar	-25°	+5°	.042	022S10-252R	022S10-252L
Single 2nd Molar	-20°	0°	.042	022S10-240R	022S10-240L
Single	0°	+5°	.042	022S10-203R	022S10-203L
Universal Single 1st Molar	0°	0°	.028	022S10-001B	022S10-001A
Universal Single	-10°	0°	.028	022S10-011B	022S10-011A

## Prescriptions

### NICKELFREE<sup>®</sup> BUCCAL TUBE



Double Buccal Tube Roth Prescription .018

Tooth	Torque	Offset	In/Out	Right	Left
<b>Upper</b>				Item #	Item #
Double 1st Molar	0°	+9°		018D10-105R	018D10-105L
Double 1st Molar	-10°	0°		018D10-110R	018D10-110L
Double 1st Molar	-10°	+9°		018D10-113R	018D10-113L
Double 1st Molar	-14°	+9°		018D10-130R	018D10-130L
Double 1st Molar	-14°	+14°		018D10-132R	018D10-132L
Universal Double	0°	0°		018D10-000A	018D10-000B
<b>Lower</b>					
Double 1st Molar	0°	+4°		018D10-202R	018D10-202L
Double 1st Molar	-20°	0°		018D10-240R	018D10-240L
Double 1st Molar	-25°	+4°		018D10-251R	018D10-251L
Universal Double	0°	0°		018D10-000B	018D10-000A



Double Buccal Tube Roth Prescription .022

Tooth	Torque	Offset	In/Out	Right	Left
<b>Upper</b>				Item #	Item #
Double 1st Molar	0°	+9°		022D10-105R	022D10-105L
Double 1st Molar	-10°	0°		022D10-110R	022D10-110L
Double 1st Molar	-10°	+9°		022D10-113R	022D10-113L
Double 1st Molar	-14°	+9°		022D10-130R	022D10-130L
Double 1st Molar	-14°	+14°		022D10-132R	022D10-132L
Universal Double	0°	0°		022D10-000A	022D10-000B
<b>Lower</b>					
Double 1st Molar	0°	+4°		022D10-202R	022D10-202L
Double 1st Molar	-20°	0°		022D10-240R	022D10-240L
Double 1st Molar	-25°	+4°		022D10-251R	022D10-251L
Universal Double	0°	0°		022D10-000B	022D10-000A

## Prescriptions

### NICKELFREE<sup>®</sup> BUCCAL TUBE



Triple Buccal Tube Roth Prescription .018


Tooth	Torque	Offset	In/Out	Right	Left
<b>Upper</b>				Item #	Item #
Triple 1st Molar	0°	+9°		018T10-105R	018T10-105L
Triple 1st Molar	-10°	+9°		018T10-113R	018T10-113L
Triple 1st Molar	-14°	+9°		018T10-130R	018T10-130L
Triple 1st Molar	-14°	+14°		018T10-132R	018T10-132L
Triple 1st Molar	-10°	+14°		018T10-114R	018T10-114L
<b>Lower</b>					



Triple Buccal Tube Roth Prescription .022

Tooth	Torque	Offset	In/Out	Right	Left
<b>Upper</b>				Item #	Item #
Triple 1st Molar	0°	+9°		022T10-105R	022T10-105L
Triple 1st Molar	-10°	+9°		022T10-113R	022T10-113L
Triple 1st Molar	-14°	+9°		022T10-130R	022T10-130L
Triple 1st Molar	-14°	+14°		022T10-132R	022T10-132L
Triple 1st Molar	-10°	+14°		022T10-114R	022T10-114L
<b>Lower</b>					

## NICKELFREE® BUCCAL TUBE



### Lip Bumper Tube Roth Prescription .018

Tooth	Torque	Offset	In/Out	Right	Left
Lower				Item #	Item #
Lip Bumper 1st Molar	0°	+4°		018T10-202R	018T10-202L
Lip Bumper 1st Molar	-12°	+4°		018T10-221R	018T10-221L
Lip Bumper 1st Molar	-25°	+4°		018T10-251R	018T10-251L

### Lip Bumper Tube Roth Prescription .022

Tooth	Torque	Offset	In/Out	Right	Left
Lower				Item #	Item #
Lip Bumper 1st Molar	0°	+4°		022T10-202R	022T10-202L
Lip Bumper 1st Molar	-12°	+4°		022T10-221R	022T10-221L
Lip Bumper 1st Molar	-25°	+4°		022T10-251R	022T10-251L



### Mini Buccal Tube Roth Prescription .018

Tooth	Torque	Offset	In/Out	Right	Left
Upper				Item #	Item #
Single - Universal	0°	0°	.023	018T10-100A	018T10-100B
Single	-10°	0°	.023	018T10-110R	018T10-110L
Lower					
Single	-12°	0°	.023	018T10-270R	018T10-270L
Single	-25°	0°	.023	018T10-280R	018T10-280L
Single - Universal	0°	0°	.023	018T10-100B	018T10-100A

### Mini Buccal Tube Roth Prescription .022

Tooth	Torque	Offset	In/Out	Right	Left
Upper				Item #	Item #
Single - Universal	0°	0°	.023	022T10-100A	022T10-100B
Single	-10°	0°	.023	022T10-110R	022T10-110L
Lower					
Single	-12°	0°	.023	022T10-270R	022T10-270L
Single	-25°	0°	.023	022T10-280R	022T10-280L
Single - Universal	0°	0°	.023	022T10-100B	022T10-100A