

STAINLESS STEEL non-convertible upper prescriptions and order info



STAINLESS STEEL NON-CONVERTIBLE UPPER BUCCAL TUBES .022.

Tooth	Torque	Offset	HDG	Right .022	Left .022	Right .018	Left .018
Single - 1st Molar	-10°	+8°	No	907.2001	907.2002	907.1001	907.1002
Single	0°	+8°	No	907.2003	907.2004	907.1003	907.1004
Single - 1st / 2nd Molar	-14°	+8°	No	907.2005	907.2006	907.1005	907.1006
Double - 1st Molar - HDG	0°	0°	.045	907.2007	907.2008	907.1007	907.1008
Double - 1st Molar	0°	+9°	No	907.2010	907.2011	907.1010	907.1011
Double - 1st Molar - HDG	0°	+9°	.045	907.2012	907.2013	907.1012	907.1013
Double - 1st Molar	-10°	0°	No	907.2014	907.2015	907.1014	907.1015
Double - 1st Molar - HDG	-10°	0°	.045	907.2016	907.2017	907.1016	907.1017
Double - 1st Molar	-10°	+9°	No	907.2018	907.2019	907.1018	907.1019
Double - 1st Molar - HDG	-10°	+9°	.045	907.2020	907.2021	907.1020	907.1021
Double - 1st Molar	-14°	+9°	No	907.2022	907.2023	907.1022	907.1023
Double - 1st Molar - HDG	-14°	+9°	.045	907.2024	907.2025	907.1024	907.1025
Double - 1st Molar	-14°	+14°	No	907.2026	907.2027	907.1026	907.1027
Double - 1st Molar - HDG	-14°	+14°	.045	907.2028	907.2029	907.1028	907.1029
Triple - 1st Molar	0°	+9°	.045	907.2030	907.2031	907.1030	907.1031
Triple - 1st Molar	-10°	+9°	.045	907.2032	907.2033	907.1032	907.1033
Triple - 1st Molar	-14°	+9°	.045	907.2034	907.2035	907.1034	907.1035
Triple - 1st Molar	-14°	+14°	.045	907.2036	907.2037	907.1036	907.1037
Triple - 1st Molar	-10°	+14°	.045	907.2038	907.2039	907.1038	907.1039

STAINLESS STEEL NON-CONVERTIBLE LOWER BUCCAL TUBES .022. & .018

Tooth	Torque	Offset	HDG	Right .022	Left .022	Right .018	Left .018
Single - 1st Molar	-25°	+5°	No	907.2040	907.2041	907.1041	907.1040
Single	-20°	0°	No	907.2042	907.2043	907.1043	907.1042
Single - 1st / 2nd Molar	0°	+5°	No	907.2044	907.2045	907.1045	907.1044
Double - 1st Molar - HDG	0°	0°	.045	907.2046	907.2047	907.1046	907.1047
Double - 1st Molar	-20°	+4°	No	907.2048	907.2049	907.1048	907.1049
Double - 1st Molar - HDG	-25°	0°	.045	907.2050	907.2051	907.1050	907.1051
Lip Bumper - 1st Molar	0°	+4°	No	907.2052	907.2053	907.1052	907.1053
Lip Bumper - 1st Molar	-12°	+4°	No	907.2054	907.2055	907.1054	907.1055
Lip Bumper - 1st Molar	-25°	+4°	No	907.2056	907.2057	907.1056	907.1057

STAINLESS STEEL NON-CONVERTIBLE UNIVERSAL BUCCAL TUBES .022

Tooth	Torque	Offset	In/out	Right/Left	.022	.018
Single - 1st Molar	0°	0°	.028	Upper Right / Lower Left	907.2058	907.1058
Single - 1st Molar	0°	0°	.028	Upper Left / Lower Right	907.2059	907.1059
Single	-10°	0°	.028	Upper Right / Lower Left	907.2060	907.1060
Single	-10°	0°	.028	Upper Left / Lower Right	907.2061	907.1061
Double	0°	0°	No	Upper Right / Lower Left	907.2009	907.1009
Double	0°	0°	No	Upper Left / Lower Right	907.2096	908.1094

STAINLESS STEEL NON-CONVERTIBLE MINI UPPER BUCCAL TUBES .022 & .018

Tooth	Torque	Offset	In/out	Right/Left	.022	.018
Single Hook	0°	0°	No	Upper Right / Lower Left	907.2062	907.1062
Single Hook	0°	0°	No	Lower Left / Upper Right	907.2063	907.1063
Single	0°	0°	No	Upper Right / Lower Left	907.2068	907.1068
Single	0°	0°	No	Lower Left / Upper Right	907.2069	907.1069
Single Hook	-10°	0°	No	Right	907.2064	907.1064
Single Hook	-10°	0°	No	Left	907.2065	907.1065
Single	-10°	0°	No	Right	907.2070	907.1070
Single	-10°	0°	No	Left	907.2071	907.1071

STAINLESS STEEL NON-CONVERTIBLE MINI LOWER BUCCAL TUBES .022 & .018

Tooth	Torque	Offset	In/out	Right .022	Left .022	Right .018	Left .018
Single Hook	-12°	0°	.023	907.2066	907.2067	907.1066	907.1067
Single	-12°	0°	.023	907.2072	907.2073	907.1072	907.1073
Single Hook	-25°	0°	.023	907.2074	907.2075	907.1074	907.1075
Single	-25°	0°	.023	907.2076	907.2077	907.1076	907.1077