

## STAINLESS STEEL BRACKETS Roth\* .022

Tooth	Torque	Offset	Item # .022
11 - Central	+12°	+5°	905.2001
12 - Lateral	+8°	+9°	905.2003
13 - Cuspid	-2°	+13°	905.2005
13 - Cuspid w/h	-2°	+13°	905.2007
14 - 1. Bicuspid	-7°	0°	905.2009
14 - 1. Bicuspid w/h	-7°	0°	905.2010
15 - 2. Bicuspid	-7°	0°	905.2009
15 - 2. Bicuspid w/h	-7°	0°	905.2010

Item # .022	Offset	Torque	Tooth
905.2002	5°	12°	Central - 21
905.2004	9°	8°	Lateral - 22
905.2006	9°	-2°	Cuspid - 23
905.2008	0°	-7°	Cuspid w/h - 23
905.2009	0°	-7°	1. Bicuspid - 24
905.2011	0°	-7°	1. Bicuspid w/h - 24
905.2009	0°	-7°	2. Bicuspid - 25
905.2011	0°	-7°	2. Bicuspid w/h - 25

Tooth	Torque	Offset	Item # .022
41 - Anterior	0°	0°	905.2014
42 - Anterior	0°	0°	905.2036
43 - Cuspid	-11°	+7°	905.2015
43 - Cuspid w/h	-11°	+7°	905.2017
44 - 1. Bicuspid	-17°	0°	905.2019
44 - 1. Bicuspid w/h	-17°	0°	905.2020
45 - 2. Bicuspid	-22°	0°	905.2022
45 - 2. Bicuspid w/h	-22°	0°	905.2023

Item # .022	Offset	Torque	Tooth
905.2014	0°	0°	Anterior - 31
905.2036	0°	0°	Anterior - 32
905.2016	+7°	-11°	Cuspid - 33
905.2018	+7°	-11°	Cuspid w/h - 33
905.2019	0°	-17°	1. Bicuspid - 34
905.2021	0°	-17°	1. Bicuspid w/h - 34
905.2022	0°	-22°	2. Bicuspid - 35
905.2024	0°	-22°	2. Bicuspid w/h - 35

## STAINLESS STEEL BRACKETS Roth\* .018

Tooth	Torque	Offset	Item # .018
11 - Central	+12°	+5°	905.1001
12 - Lateral	+8°	+9°	905.1003
13 - Cuspid	-2°	+13°	905.1005
13 - Cuspid w/h	-2°	+13°	905.1007
14 - 1. Bicuspid	-7°	0°	905.1009
14 - 1. Bicuspid w/h	-7°	0°	905.1010
15 - 2. Bicuspid	-7°	0°	905.1009
15 - 2. Bicuspid w/h	-7°	0°	905.1010

Item # .018	Offset	Torque	Tooth
905.1002	5°	12°	Central - 21
905.1004	9°	8°	Lateral - 22
905.1006	9°	-2°	Cuspid - 23
905.1008	0°	-7°	Cuspid w/h - 23
905.1009	0°	-7°	1. Bicuspid - 24
905.1011	0°	-7°	1. Bicuspid w/h - 24
905.1009	0°	-7°	2. Bicuspid - 25
905.1011	0°	-7°	2. Bicuspid w/h - 25

Tooth	Torque	Offset	Item # .018
41 - Anterior	0°	0°	905.1014
42 - Anterior	0°	0°	905.1014
43 - Cuspid	-11°	+7°	905.1015
43 - Cuspid w/h	-11°	+7°	905.1017
44 - 1. Bicuspid	-17°	0°	905.1019
44 - 1. Bicuspid w/h	-17°	0°	905.1020
45 - 2. Bicuspid	-22°	0°	905.1022
45 - 2. Bicuspid w/h	-22°	0°	905.1023

Item # .018	Offset	Torque	Tooth
905.1014	0°	0°	Anterior - 31
905.1036	0°	0°	Anterior - 32
905.1016	+7°	-11°	Cuspid - 33
905.1018	+7°	-11°	Cuspid w/h - 33
905.1019	0°	-17°	1. Bicuspid - 34
905.1021	0°	-17°	1. Bicuspid w/h - 34
905.1022	0°	-22°	2. Bicuspid - 35
905.1024	0°	-22°	2. Bicuspid w/h - 35